

CytoScientific Centrifuge Tubes

CytoScientific





Certificate of Registration

This is to certify that

CYTO SCIENTIFIC

2ND FLOOR, SCF NO 38 SINGLE SEAT, GURU NANAK SWEETS, PHASE 10,
MOHALLI, SAS NAGAR, PUNJAB, 160055, INDIA

has been independently assessed by QRO
and is compliant with the requirement of:

ISO 13485:2016

Quality Management System for Medical Devices

For the following scope of activities:

MANUFACTURING OF UNIVERSAL PIPETTE TIPS, MICROCENTRIFUGE TUBES, CENTRIFUGE TUBES, PCR TUBES, 8 PCR TUBES STRIP, PCR PLATE, DEEP WELL PLATE, LABORATORY CENTRIFUGE, CRYOVIALS, CRYO BOX, NITRILE GLOVES, VINYL GLOVES, LATEX EXAMINATION GLOVES, FACE MASK, SURGEON CAP, SHOE COVER, DISPOSABLE MEDICAL GOWN, DISPOSABLE MEDICAL APRON, PERSONAL PROTECTIVE CLOTHING (PPE) KIT, SEXUAL ASSAULT EVIDENCE COLLECTION KIT (SAECK), FILTER PAPER, DNA STORAGE CARD, BLOOD COLLECTION CARDS, STERILE FOAM TIPPED APPLICATOR, COTTON-TIPPED SWAB

Date of Certification: 20th October 2022

2nd Surveillance Audit Due: 19th October 2024

1st Surveillance Audit Due: 19th October 2023

Certificate Expiry: 19th October 2025

Certificate Number: 305022102045MD



Chauhan

Head of Certification

EGAC is a member of International Accreditation forum (IAF) for QMS, EMS, OHSM, FSMS and EnMS

Validity of this certificate is subject to annual surveillance audits to be done successfully on or before 365 days from date of the audit. (In case surveillance audit is not allowed to be conducted: this certificate shall be suspended / withdrawn).

The Validity of this certificate can be verified at www.grocert.org

This certificate of registration remains the property of QRO Certification LLP, and shall be returned immediately upon request.

India Office : QRO Certification LLP

142, IInd Floor, Avtar Enclave, Near Paschim Vihar West Metro Station, Delhi-110063, (INDIA)

Website : www.grocert.org, E-mail : info@grocert.org

2/4/22, 2:52 PM

Print : Udyam Registration Certificate



भारत सरकार
Government of India
सूक्ष्म, लघु एवं मध्यम उद्यम मंत्रालय
Ministry of Micro, Small and Medium Enterprises



UDYAM REGISTRATION CERTIFICATE



UDYAM REGISTRATION NUMBER	UDYAM-PB-20-0015551				
NAME OF ENTERPRISE	CYTO SCIENTIFIC				
TYPE OF ENTERPRISE *	MICRO				
MAJOR ACTIVITY	MANUFACTURING				
SOCIAL CATEGORY OF ENTREPRENEUR	GENERAL				
NAME OF UNIT(S)	S.No.	Udyog Aadhaar Memorandum	Unit(s) Name		
	1	PR20D0013179	Cyto Scientific		
OFFICIAL ADDRESS OF ENTERPRISE	Flat/Door/Block No.	SCF 38	Name of Premises/ Building	II Floor	
	Village/Town	Mohali	Block	Phase 10	
	Road/Street/Lane	Near Guru Nanak Sweets	City	Mohali	
	State	PUNJAB	District	SAS NAGAR , Pin 160062	
	Mobile	955784826	Email:	cytoscientific@gmail.com	
DATE OF INCORPORATION / REGISTRATION OF ENTERPRISE	05/09/2019				
DATE OF COMMENCEMENT OF PRODUCTION/BUSINESS	12/04/2021				
NATIONAL INDUSTRY CLASSIFICATION CODE(S)	S.No.	NIC 2 Digit	NIC 4 Digit	NIC 5 Digit	Activity
	1	21 - Manufacture of pharmaceuticals, medicinal chemical and botanical products	2100 - Manufacture of pharmaceuticals, medicinal chemical and botanical products	21006 - Manufacture of medical impregnated wadding, gauze, bandages, dressings, surgical gut string etc.	Manufacturing
	2	22 - Manufacture of rubber and plastics products	2220 - Manufacture of plastics products	22209 - Manufacture of other plastics products n.e.e	Manufacturing
	3	28 - Manufacture of machinery and equipment n.e.e.	2829 - Manufacture of other special-purpose machinery	28299 - Manufacture of other special-purpose machinery n.e.e.	Manufacturing
4	74 - Other professional, scientific and technical activities n.e.e.	7490 - Other professional, scientific and technical activities n.e.e.	74901 - Business brokerage activities	Services	
DATE OF UDYAM REGISTRATION	18/09/2021				

* In case of graduation (upward/reverse) of status of an enterprise, the benefit of the Government Schemes will be availed as per the provisions of Notification No. S.O. 2119(E) dated 28.06.2020 issued by the Mo/MSME.

Disclaimer: This is computer generated statement, no signature required. Printed from <https://udyamregistration.gov.in> & Date of printing: 04/02/2022

For any assistance, you may contact:

1. District Industries Centre: SAS NAGAR (PUNJAB)

2. MSME-DI: LUDHIANA (PUNJAB)

<https://udyamregistration.gov.in/PrintApplication.aspx?fdm=6rRslYu+cZvSbLLprq2UA==>

BE A
CHAMPION
with the
Ministry of


1/4

CytoScientific Centrifuge Tubes are designed to meet the rigorous demands of your laboratory applications. Available in various sizes and compatible with most centrifuges, these tubes ensure safety, durability, and consistent performance.

Volume - 15ml




Volume - 50ml



Conical Bottom

Key Features Of Centrifuge Tubes



Legible Mark

Grduated Markes



Leak-proof Screw Cap

Conical Bottom

- Made of Medical grade virgin polypropylene
- Color- Transparent or clear
- Conical Bottom
- Autoclavable
- Graduated
- Sterile
- Large Surface area for labelling
- Tubes are thermal & freezing stable, chemical resistant, non pyrogenic, non-cytotoxic.

Order Your Centrifuge Tubes Today

Sr. No.	Volume	Pack Size	Catalog No.
1.	15ml	50 Tubes per pack	CS-CT-15ml-S
2.	50ml	25 Tubes per pack	CS-CT-50ml-S

For more information, contact us or visit our website.

Phone No. - [+91-8130681839](tel:+91-8130681839), [+91-7696220884](tel:+91-7696220884)

Email - cytoscientific@gmail.com

Website - <https://cytoscientific.com/>

