

# CytoScientific Cryovials

CytoScientific





# Certificate of Registration

This is to certify that

## CYTO SCIENTIFIC

2ND FLOOR, SCF NO 38 SINGLE SEAT, GURU NANAK SWEETS, PHASE 10,  
MOHALI, SAS NAGAR, PUNJAB, 160055, INDIA

has been independently assessed by QRO  
and is compliant with the requirement of:

**ISO 13485:2016**

**Quality Management System for Medical Devices**

For the following scope of activities:

MANUFACTURING OF UNIVERSAL PIPETTE TIPS, MICROCENTRIFUGE TUBES, CENTRIFUGE TUBES, PCR TUBES, 8 PCR TUBES STRIP, PCR PLATE, DEEP WELL PLATE, LABORATORY CENTRIFUGE, CRYOVIALS, CRYO BOX, NITRILE GLOVES, VINYL GLOVES, LATEX EXAMINATION GLOVES, FACE MASK, SURGEON CAP, SHOE COVER, DISPOSABLE MEDICAL GOWN, DISPOSABLE MEDICAL APRON, PERSONAL PROTECTIVE CLOTHING (PPE) KIT, SEXUAL ASSAULT EVIDENCE COLLECTION KIT (SAECK), FILTER PAPER, DNA STORAGE CARD, BLOOD COLLECTION CARDS, STERILE FOAM TIPPED APPLICATOR, COTTON-TIPPED SWAB

Date of Certification: 20th October 2022

2<sup>nd</sup> Surveillance Audit Due: 19th October 2024

1<sup>st</sup> Surveillance Audit Due: 19th October 2023

Certificate Expiry: 19th October 2025

**Certificate Number: 305022102045MD**



*Chauran*  
Head of Certification

EGAC is a member of International Accreditation forum (IAF) for QMS, EMS, OHSM, FSMS and EnMS

Validity of this certificate is subject to annual surveillance audits to be done successfully on or before 365 days from date of the audit.  
(In case surveillance audit is not allowed to be conducted: this certificate shall be suspended / withdrawn).

The Validity of this certificate can be verified at [www.qrocert.org](http://www.qrocert.org)

This certificate of registration remains the property of QRO Certification LLP, and shall be returned immediately upon request.

India Office : QRO Certification LLP

142, IInd Floor, Astar Enclave, Near Paschim Vihar West Metro Station, Delhi-110063, (INDIA)

Website : [www.qrocert.org](http://www.qrocert.org), E-mail : [info@qrocert.org](mailto:info@qrocert.org)

2/4/22, 2:52 PM

Print : Udyam Registration Certificate



भारत सरकार  
Government of India  
सूक्ष्म, लघु एवं मध्यम उद्यम मंत्रालय  
Ministry of Micro, Small and Medium Enterprises



## UDYAM REGISTRATION CERTIFICATE



UDYAM REGISTRATION NUMBER	UDYAM-PB-20-0015551			
NAME OF ENTERPRISE	CYTO SCIENTIFIC			
TYPE OF ENTERPRISE *	MICRO			
MAJOR ACTIVITY	MANUFACTURING			
SOCIAL CATEGORY OF ENTREPRENEUR	GENERAL			
NAME OF UNIT(S)	S.No. Udyog Aadhaar Memorandum 1 PB200013179	Unit(s) Name Cyto Scientific		
OFFICIAL ADDRESS OF ENTERPRISE	Flat/Door/Block No. SCF 38 Village/Town Mohali Road/Street/Lane Near Guru Nanak Sweets State PUNJAB Mobile 9555784836	Name of Premises/ Building Block City Mohali District SAS NAGAR, Pin 160062 Email: cytoscientific@gmail.com	II Floor Phase 10	
DATE OF INCORPORATION / REGISTRATION OF ENTERPRISE	05/09/2019			
DATE OF COMMENCEMENT OF PRODUCTION/BUSINESS	12/04/2021			
NATIONAL INDUSTRY CLASSIFICATION CODE(S)	S.No. NIC 2 Digit 1 21 - Manufacture of pharmaceuticals, medicinal chemical and botanical products 2 22 - Manufacture of rubber and plastics products 3 28 - Manufacture of machinery and equipment n.e.c. 4 74 - Other professional, scientific and technical activities n.e.c.	NIC 4 Digit 2100 - Manufacture of pharmaceuticals, medicinal chemical and botanical products 2220 - Manufacture of plastics products 2829 - Manufacture of other special-purpose machinery 7490 - Other professional, scientific and technical activities n.e.c.	NIC 5 Digit 21006 - Manufacture of medical impregnated wadding, gauze, bandages, dressings, surgical gut string etc. 22299 - Manufacture of other plastics products n.e.c. 28299 - Manufacture of other special-purpose machinery n.e.c. 74901 - Business brokerage activities n.e.c.	Activity Manufacturing Manufacturing Manufacturing Services
DATE OF UDYAM REGISTRATION	18/09/2021			

\* In case of graduation (upward/reverse) of status of an enterprise, the benefit of the Government Schemes will be availed as per the provisions of Notification No. S.O. 2119(E) dated 26.06.2020 issued by the Mo/MSME.

Disclaimer: This is computer generated statement, no signature required. Printed from <https://udyamregistration.gov.in> & Date of printing: 04/02/2022

For any assistance, you may contact:

1. District Industries Centre: SAS NAGAR (PUNJAB)

2. MSME-DI: LUDHIANA (PUNJAB)

**BE A  
CHAMPION  
with the  
Ministry of**


<https://udyamregistration.gov.in/PrintApplication.aspx?udm=6rRslYU+cZvSbLLprq2UA==>

1/4

Cryovials are essential tools for laboratories that require the secure, long-term storage of biological samples. Whether you're working in biobanking, research, or clinical studies, our cryovials ensure your samples are preserved safely at ultra-low temperatures.

## External Thread Cryovials

- Color-coded polypropylene caps for easy identification.
- Silicone washer provides a tight, leak-proof seal.
- Clear black graduations for easy volume measurement.
- Stable, self-standing base for convenience.
- Self-locking skirt for added security.



## Internal Thread Cryovials

- Color-coded polypropylene caps for easy identification.
- Silicone washers or rubber O-rings ensure a leak-proof seal.
- Easy-to-read black graduations for accurate measurements.
- Self-standing base and self-locking skirt for convenience.

# Ordering Information



- Made from high-quality polypropylene.
- Withstands temperatures as low as  $-196^{\circ}\text{C}$  (ideal for liquid nitrogen).
- Leak-proof design with secure screw cap and silicone seal.
- Self-standing base for stability during storage and transport.
- Easy one-handed handling with a wide mouth for filling/retrieval.
- Frosted area for easy and permanent labeling.
- RNase/DNase-free and non-pyrogenic to minimize contamination risk.

Volume	Pack Size	Catalog No.
1.8ml/ 2.0 ml	100/500 Tubes per pack	CS-CV-2.0ml-S
5.0 ml	100 Tubes per pack	CS-CV-5.0ml-S

## Cryovials Storage Boxes



- Made from virgin polypropylene with 99.99% purity, suitable for use in IP grade applications.
- With a lid featuring alphanumeric labeling, it can be used from  $-90^{\circ}\text{C}$  to  $121^{\circ}\text{C}$ .
- Autoclavable
- Available in various sizes and configurations to fit different numbers of tubes (81, 96 or 100).

Item	No. of Holes	Catalog No.
Cryobox - 81	81	CS-CB-81
Cryobox - 100	100	CS-CB-100

# Order Your Cryovials Today!

Ensure the safety and integrity of your valuable samples with our high-quality cryovials. Contact us today for more information or to place an order.

Phone No. - +91-8130681839, +91-7696220884

Email - [cytoscientific@gmail.com](mailto:cytoscientific@gmail.com)

Website - <https://cytoscientific.com/>