


# CytoScientific Deep well plates




Deep well plates are laboratory tools used for storing and handling large volumes of liquid samples. Made from durable plastics like polypropylene, they are resistant to chemicals and can be used in various applications such as high-throughput screening, sample storage, and mixing.



- Made of polypropylene
- Dnase/Rnase free
- Compatible with Automated system
- Sealing film
- Autoclavable
- Suitable for multi channel pipettes

Material	Polypropylene
Features	Dnase/Rnase free ,Compatible with Automated system ,Sealing film ,Autoclavable ,Suitable for multi channel pipettes
Colour	White/Clear
No Of wells	96
Type of Bottom	U well / V well
Number of piece in one box	100



# Order Information

Item	Volume	Pack Size	Catalog No.
96-Deep-Well Plate Square Well Round Conical Bottom (Elution Plate)	200ul	50	CS-96-0.5-RC
96-Deep-Well Plate Square Well Round Conical Bottom (Elution Plate)	1.6ml	50	CS-96-1.6-RC
96-Deep-Well Plate Square Well Round Conical Bottom (Elution Plate)	2.2ml	50	CS-96-2.2-SC
96-Deep-Well Plate Square Well Round U-Bottom (Elution Plate)	2.2ml	50	CS-96-2.2-SU

## Contact Us

For inquiries or to place an order, reach out to us today!

Email: [cytoscientific@gmail.com](mailto:cytoscientific@gmail.com)

Phone: [+91-8130681839](tel:+91-8130681839), [+91-7696220884](tel:+91-7696220884)

Website: [www.cytoscientific.com](http://www.cytoscientific.com)