


CytoScientific Microtiter Plate

CytoScientific





Certificate of Registration

This is to certify that

CYTO SCIENTIFIC

2ND FLOOR, SCF NO 38 SINGLE SEAT, GURU NANAK SWEETS, PHASE 10,
MOHALI, SAS NAGAR, PUNJAB, 160055, INDIA

has been independently assessed by QRO
and is compliant with the requirement of:

ISO 13485:2016

Quality Management System for Medical Devices

For the following scope of activities:

MANUFACTURING OF UNIVERSAL PIPETTE TIPS, MICROCENTRIFUGE TUBES, CENTRIFUGE TUBES, PCR TUBES, 8 PCR TUBES STRIP, PCR PLATE, DEEP WELL PLATE, LABORATORY CENTRIFUGE, CRYOVIALS, CRYO BOX, NITRILE GLOVES, VINYL GLOVES, LATEX EXAMINATION GLOVES, FACE MASK, SURGEON CAP, SHOE COVER, DISPOSABLE MEDICAL GOWN, DISPOSABLE MEDICAL APRON, PERSONAL PROTECTIVE CLOTHING (PPE) KIT, SEXUAL ASSAULT EVIDENCE COLLECTION KIT (SAECK), FILTER PAPER, DNA STORAGE CARD, BLOOD COLLECTION CARDS, STERILE FOAM TIPPED APPLICATOR, COTTON-TIPPED SWAB

Date of Certification: 20th October 2022

2nd Surveillance Audit Due: 19th October 2024

1st Surveillance Audit Due: 19th October 2023

Certificate Expiry: 19th October 2025

Certificate Number: 305022102045MD



Chaurhary
Head of Certification

EGAC is a member of International Accreditation forum (IAF) for QMS, EMS, OHSM, FSMS and EnMS

Validity of this certificate is subject to annual surveillance audits to be done successfully on or before 365 days from date of the audit.
(In case surveillance audit is not allowed to be conducted: this certificate shall be suspended / withdrawn).

The Validity of this certificate can be verified at www.qrocert.org

This certificate of registration remains the property of QRO Certification LLP, and shall be returned immediately upon request.

India Office : QRO Certification LLP

142, IInd Floor, Astar Enclave, Near Paschim Vihar West Metro Station, Delhi-110063, (INDIA)

Website : www.qrocert.org, E-mail : info@qrocert.org

2/4/22, 2:52 PM

Print : Udyam Registration Certificate



भारत सरकार
Government of India
सूक्ष्म, लघु एवं मध्यम उद्यम मंत्रालय
Ministry of Micro, Small and Medium Enterprises



UDYAM REGISTRATION CERTIFICATE



UDYAM REGISTRATION NUMBER	UDYAM-PB-20-0015551			
NAME OF ENTERPRISE	CYTO SCIENTIFIC			
TYPE OF ENTERPRISE *	MICRO			
MAJOR ACTIVITY	MANUFACTURING			
SOCIAL CATEGORY OF ENTREPRENEUR	GENERAL			
NAME OF UNIT(S)	S.No. Udyog Aadhaar Memorandum 1 PB20D0013179	Unit(s) Name Cyto Scientific		
OFFICIAL ADDRESS OF ENTERPRISE	Flat/Door/Block No. SCF 38 Village/Town Mohali Road/Street/Lane Near Guru Nanak Sweets State PUNJAB Mobile 9555784836	Name of Premises/ Building Block City Mohali District SAS NAGAR, Pin 160062 Email: cytoscientific@gmail.com	II Floor Phase 10	
DATE OF INCORPORATION / REGISTRATION OF ENTERPRISE	05/09/2019			
DATE OF COMMENCEMENT OF PRODUCTION/BUSINESS	12/04/2021			
NATIONAL INDUSTRY CLASSIFICATION CODE(S)	S.No. NIC 2 Digit 1 21 - Manufacture of pharmaceuticals, medicinal chemical and botanical products 2 22 - Manufacture of rubber and plastics products 3 28 - Manufacture of machinery and equipment n.e.c. 4 74 - Other professional, scientific and technical activities n.e.c.	NIC 4 Digit 2100 - Manufacture of pharmaceuticals, medicinal chemical and botanical products 2220 - Manufacture of plastics products 2829 - Manufacture of other special-purpose machinery 7490 - Other professional, scientific and technical activities n.e.c.	NIC 5 Digit 21006 - Manufacture of medical impregnated wadding, gauze, bandages, dressings, surgical gut string etc. 22299 - Manufacture of other plastics products n.e.c. 28299 - Manufacture of other special-purpose machinery n.e.c. 74901 - Business brokerage activities n.e.c.	Activity Manufacturing Manufacturing Manufacturing Services
DATE OF UDYAM REGISTRATION	18/09/2021			

* In case of graduation (upward/reverse) of status of an enterprise, the benefit of the Government Schemes will be availed as per the provisions of Notification No. S.O. 2119(E) dated 26.06.2020 issued by the Mo/MSME.

Disclaimer: This is computer generated statement, no signature required. Printed from <https://udyamregistration.gov.in> & Date of printing: 04/02/2022

For any assistance, you may contact:

1. District Industries Centre: SAS NAGAR (PUNJAB)




2. MSME-DI: LUDHIANA (PUNJAB)

<https://udyamregistration.gov.in/PrintApplication.aspx?udm=6rRslYU+cZvSbLLprq2UA==>


**BE A
CHAMPION
with the
Ministry of**

Microtiter plates are flat, multi-well plates used in laboratories for various assays and experiments. Microtiter plates, also known as microplates or multiwell plates, are compact, flat plates feature an array of tiny wells, functioning like miniature test tubes for high-throughput applications.

- Made of polystyrene and available in sterile or non – sterile options
- Available in various well counts including 6, 24, 96 or 384 wells.
- Different color options – clear, white or black, depending on their intended use.
- Designed for easy stacking and storage.
- Available with lids to prevent contamination
- Different well surface treatments are available for cell growth or non-stick surfaces.
- Barcode option also available for easy identification and tracking.
- Compatible with automatic pipetting systems, plate readers, and other laboratory instruments



Material	Polypropylene
Color	Clear, White or Black
Plate size	6, 24, 96 or 384 wells
Bottom	U bottom,V bottom,Flat bottom for choice
RNase- and DNase-free	Yes



Order Your Microtiter plate Today

Item	Pack Size	Catalog No.
Microtitration Plates, Flat Bottom	100	CS-MTP-96F-S
Microtitration Plates, Round Bottom	100	CS-MTP-96U-S
Microtitration Plates with Lids, Flat Bottom	100	CS-MTP-96F-L-S
Microtitration Plates with Lids, Round Bottom	100	CS-MTP-96U-L-S

For more information, contact us or visit our website.

Phone No. - [+91-8130681839](tel:+91-8130681839), [+91-7696220884](tel:+91-7696220884)

Email - cytoscientific@gmail.com

Website - <https://cytoscientific.com/>