


CytoScientific PCR Plates

CytoScientific





Certificate of Registration

This is to certify that

CYTO SCIENTIFIC

2ND FLOOR, SCF NO 38 SINGLE SEAT, GURU NANAK SWEETS, PHASE 10,
MOHALI, SAS NAGAR, PUNJAB, 160055, INDIA

has been independently assessed by QRO
and is compliant with the requirement of:

ISO 13485:2016

Quality Management System for Medical Devices

For the following scope of activities:

MANUFACTURING OF UNIVERSAL PIPETTE TIPS, MICROCENTRIFUGE TUBES, CENTRIFUGE TUBES, PCR TUBES, 8 PCR TUBES STRIP, PCR PLATE, DEEP WELL PLATE, LABORATORY CENTRIFUGE, CRYOVIALS, CRYO BOX, NITRILE GLOVES, VINYL GLOVES, LATEX EXAMINATION GLOVES, FACE MASK, SURGEON CAP, SHOE COVER, DISPOSABLE MEDICAL GOWN, DISPOSABLE MEDICAL APRON, PERSONAL PROTECTIVE CLOTHING (PPE) KIT, SEXUAL ASSAULT EVIDENCE COLLECTION KIT (SAECK), FILTER PAPER, DNA STORAGE CARD, BLOOD COLLECTION CARDS, STERILE FOAM TIPPED APPLICATOR, COTTON-TIPPED SWAB

Date of Certification: 20th October 2022

2nd Surveillance Audit Due: 19th October 2024

1st Surveillance Audit Due: 19th October 2023

Certificate Expiry: 19th October 2025

Certificate Number: 305022102045MD



Chauhan

Head of Certification

EGAC is a member of International Accreditation forum (IAF) for QMS, EMS, OHSM, FSMS and EnMS

Validity of this certificate is subject to annual surveillance audits to be done successfully on or before 365 days from date of the audit. (In case surveillance audit is not allowed to be conducted: this certificate shall be suspended / withdrawn).

The Validity of this certificate can be verified at www.qrocert.org

This certificate of registration remains the property of QRO Certification LLP, and shall be returned immediately upon request.

India Office : QRO Certification LLP

142, 11nd Floor, Avtar Enclave, Near Paschim Vihar West Metro Station, Delhi-110063, (INDIA)

Website : www.qrocert.org, E-mail : info@qrocert.org

2/4/22, 2:52 PM

Print : Udyam Registration Certificate



भारत सरकार
Government of India
सूक्ष्म, लघु एवं मध्यम उद्यम मंत्रालय
Ministry of Micro, Small and Medium Enterprises



UDYAM REGISTRATION CERTIFICATE



UDYAM REGISTRATION NUMBER	UDYAM-PB-20-0015551				
NAME OF ENTERPRISE	CYTO SCIENTIFIC				
TYPE OF ENTERPRISE*	MICRO				
MAJOR ACTIVITY	MANUFACTURING				
SOCIAL CATEGORY OF ENTREPRENEUR	GENERAL				
NAME OF UNIT(S)	S.No.	Udyog Aadhaar Memorandum	Unit(s) Name		
	1	PR2000013179	Cyto Scientific		
OFFICAL ADDRESS OF ENTERPRISE	Flat/Door/Block No.	SCF 38	Name of Premises/ Building		
	Village/Town	Mohali	Block		
	Road/Street/Lane	Near Guru Nanak Sweets	City		
	State	PUNJAB	District		
	Mobile	9555784836	Email:		
			cytoscientific@gmail.com		
DATE OF INCORPORATION / REGISTRATION OF ENTERPRISE	05/09/2019				
DATE OF COMMENCEMENT OF PRODUCTION/BUSINESS	12/04/2021				
NATIONAL INDUSTRY CLASSIFICATION CODE(S)	S.No.	NIC 2 Digit	NIC 4 Digit	NIC 5 Digit	Activity
	1	21 - Manufacture of pharmaceuticals, medicinal chemical and botanical products	2100 - Manufacture of pharmaceuticals, medicinal chemical and botanical products	21006 - Manufacture of medicinal impregnated wadding, gauze, bandages, dressings, surgical gut string etc.	Manufacturing
	2	22 - Manufacture of rubber and plastics products	2220 - Manufacture of plastics products	22209 - Manufacture of other plastics products n.e.c.	Manufacturing
	3	28 - Manufacture of machinery and equipment n.e.c.	2829 - Manufacture of other special-purpose machinery	28299 - Manufacture of other special-purpose machinery n.e.c.	Manufacturing
	4	74 - Other professional, scientific and technical activities n.e.c.	7490 - Other professional, scientific and technical activities n.e.c.	74901 - Business brokerage activities	Services
DATE OF UDYAM REGISTRATION	18/09/2021				

* In case of graduation (upward/reverse) of status of an enterprise, the benefit of the Government Schemes will be availed as per the provisions of Notification No. S.O. 2119(E) dated 26.06.2020 issued by the MoU MSME.

Disclaimer: This is computer generated statement, no signature required. Printed from <https://udyamregistration.gov.in> & Date of printing:- 04/02/2022

For any assistance, you may contact:

- District Industries Centre: SAS NAGAR (PUNJAB)
- MSME-DI: LUDHIANA (PUNJAB)

BE A
CHAMPION
with the
Ministry of

<https://udyamregistration.gov.in/PrintApplication.aspx?Udm=6rRsllyU+cZvSbLLprq2UA==>

1/4

PCR plates, also known as PCR microplates or PCR trays, are used for conducting Polymerase Chain Reaction (PCR) experiments. These plates are designed to hold PCR reaction mixtures and allow for efficient amplification of DNA or RNA samples.




- Ultra-pure polypropylene for low contamination.
- Thin walls ensure even heat transfer.
- Clear for easy sample visibility.
- Available in 96-well or 384-well formats.
- Alphanumeric labeling for identification.
- Fits most thermal cyclers.
- Sterile, RNase- and DNase-free options.
- Skirted, semi-skirted, and non-skirted versions available.

High-contrast, black alphanumeric indexing.



Raised wells ensure optimal sealing to prevent evaporation.



Skirt options




Full-Skirted




Non-Skirted


Ultra-clear wells for easy sample visualization.



Standard Well
0.2 mL



Low-Profile
0.1 mL



Semi-Skirted

1118642
1118641
1118640
1118639
1118638

Custom bar-coding available.

Order Information

Item	Volume/Specification	Pack Size	Catalog No.
96-Well PCR PlateNo Skirt	0.1ml, Transparent or Clear	10	CS-96-010-PCR-NS-C
96-Well PCR PlateNo Skirt	0.1ml, White	10	CS-96-010-PCR-NS-W
96-Well PCR PlateNo Skirt	0.2ml, Transparent or Clear	10	CS-96-010-PCR-NS-C
96-Well PCR PlateNo Skirt	0.2ml, White	10	CS-96-010-PCR-NS-W
96-Well PCR PlateSemi Skirt	0.2ml, Transparent or Clear	10	CS-96-020-PCR-SS-C
96-Well PCR PlateSemi Skirt	0.2ml, Black	10	CS-96-020-PCR-SS-B
96-Well PCR PlateSemi Skirt	0.2ml, White	10	CS-96-020-PCR-SS-W
96-Well PCR PlateFull Skirt	0.1ml, Transparent or Clear	10	CS-96-020-PCR-FS-C

Contact Us

For inquiries or to place an order, reach out to us today!

Email: cytoscientific@gmail.com

Phone: [+91-8130681839](tel:+91-8130681839), [+91-7696220884](tel:+91-7696220884)

Website: www.cytoscientific.com