

# CYTOSCIENTIFIC

## LAB PLASTICWARE



# ABOUT US

We are a leading manufacturer of high-quality laboratory plasticware, dedicated to providing researchers, labs, and educational institutions with the reliable and innovative equipment they need to push the boundaries of knowledge.

## A Legacy of Excellence

Our story began with a commitment to quality and a vision for the future of scientific research. Since 2019, we've been at the forefront of laboratory plasticware development, continuously refining our processes and materials to deliver products that are:


- **Durable and dependable:** We understand the critical role our products play in your experiments. That's why we use only the highest-grade materials and rigorous manufacturing practices to ensure our plasticware withstands the demands of even the most challenging applications.
- **Functional and innovative:** We stay abreast of the evolving needs of scientific research, constantly developing new and improved plasticware solutions. Our team of experts works tirelessly to create products that enhance efficiency, accuracy, and safety in the lab.
- **Cost-effective:** We believe that groundbreaking science shouldn't be hindered by budgetary constraints. We offer competitive pricing without compromising on quality, ensuring you have the tools you need to achieve your research goals.

## Beyond Products: A Commitment to Partnership

At Cytoscientific, we view ourselves as more than just a manufacturer. We are a trusted partner in your scientific journey. Our knowledgeable team is available to assist you in choosing the right plasticware for your specific needs, and we offer ongoing support to ensure you get the most out of our products.

Join us in Shaping the Future of Science

# Table Of Contents



**01**

## **Pipette Tips**

- Universal Tips
- Aerosol Filter Barrier Tips

**02**

## **Tubes**

- Centrifuge Tubes
- Microcentrifuge Tubes
- Storage/Transparent Tubes

**03**

## **PCR Consumables**

- PCR Tubes/Strips
- PCR Plates
- Plate Sealers

**04**

## **Assay/Storage Plats**

- Deep well Plates
- Microtiter Plates

**05**

## **Cryogenic Vials**

- Cryogenic Vial

**06**

## **Labware**

- Sample Container
- Bottles
- Microcentrifuge Tube Racks
- General Labware

**07**

## **Sexual Assault Evidence Collection Kit**

# 1 UNIVERSAL PIPETTE TIPS

- Made of Medical grade virgin polypropylene
- Universal Fit
- Autoclavable
- Graduated
- Sterile/Non-Sterile options available



Sr. No.	Volume	Pack Size	Catalog No.
1.	0.5-10ul	Bulk (1000 Tips per pack)Color - Transparent	CS-UT-10-C
2.	2-200ul	Bulk (1000 Tips per pack)Color - Transparent, Yellow	CS-UT-200-C
3.	100-1000ul	Bulk (1000 Tips per pack)Color - Transparent,Blue	CS-UT-1000-C

# 2 Aerosol Filter Barrier Tips

- Made of medical grade virgin polypropylene
- Color- Transparent or clear
- Universal Fit
- Autoclavable
- Graduated
- Sterile
- RNase DNase and pyrogen free option also available



Sr. No.	Volume	Pack Size	Catalog No.
1.	0.5 -10ul	Bulk ( 1000Tips)Rack (96Tips)	CS-FB-10-SCS-FB-10-R-S
2.	2 -20ul	Bulk ( 1000Tips)Rack (96Tips)	CS-FB-20-SCS-FB-20-R-S
3.	10 -100ul	Bulk ( 1000Tips)Rack (96Tips)	CS-FB-100-SCS-FB-100-R-S
4.	20-200ul	Bulk ( 1000Tips)Rack (96Tips)	CS-FB-200-SCS-FB-200-R-S
5.	100-1000ul	Bulk ( 1000Tips)Rack (96Tips)	CS-FB-1000-SCS-FB-1000-R-S

# 3 Centrifuge Tubes

- Made of Medical grade virgin polypropylene
- Color- Transparent or clear
- Conical Bottom
- Autoclavable
- Graduated
- Sterile
- Large Surface area for labelling
- Tubes are thermal & freezing stable, chemical resistant, non pyrogenic, non-cytotoxic.



Sr. No.	Volume	Pack Size	Catalog No.
1.	15ml	50 Tubes per pack	CS-CT-15ml-S
2.	50ml	25 Tubes per pack	CS-CT-50ml-S

# 4 Microcentrifuge Tubes

- Made of Medical grade virgin polypropylene
- Color- Transparent or clear
- Autoclavable
- Graduated
- Sterile
- Frosted surface area for labelling
- 0.5 ml and 1.5 ml available in conical bottom, 2 ml available in conical and round bottom options
- Tubes are thermal & freezing stable, chemical resistant, non pyrogenic, non-cytotoxic.



Sr. No.	Volume	Pack Size	Catalog No.
1.	0.5ml	1000 Tubes per pack	CS-MT-060-C-S
2.	1.5ml	500 Tubes per pack	CS-MT-150-C-S
3.	2.0ml	500 Tubes per pack	CS-MT-200-C-S

## 5 Storage/Transport Tubes

- Made from autoclavable medical-grade polypropylene.
- Leak-proof with screw-top lids and secure seal.
- Sizes available: 1ml, 1.5ml, 2.0ml, 5.0ml, and 10.0ml.
- Available in sterile and non-sterile options.
- Easy labeling with clear walls and labeling areas.
- Cost-effective for sample organization.



Sr. No.	Volume	Pack Size	Catalog No.
1.	1.0ml	100 Tubes per pack	CS-SV-1.0ml
2.	1.5ml	100 Tubes per pack	CS-SV-1.5ml
3.	2.0ml	100 Tubes per pack	CS-SV-2.0ml
4.	5.0ml	100 Tubes per pack	CS-SV-5.0ml
5.	10.0ml	100 Tubes per pack	CS-SV-10.0ml

## 6 PCR Tubes Strip


- Ultra-pure, medical-grade polypropylene.
- Thin walls for efficient heat transfer and uniform PCR.
- Clear for easy sample visualization.
- Tight-fitting caps to reduce sample loss.
- Available in sterile formats: individual tubes or 4/8-tube strips.
- Volume options: 0.1mL or 0.2mL.
- Flat or domed caps available.
- RNase- and DNase-free options for sensitive work.



Item	Volume	Pack Size	Catalog No.
Individual PCR tube	0.2 ml	1000 Tubes per pack	CS-02-C-PCR
8 PCR tubes strip, Caps Separate	0.2 ml	125 Strips per pack	CS-0208-FCC-PCR
8 PCR tubes strip, Caps Attached	0.2 ml	125 Strips per pack	CS-0208-FCCA-PCR

# 7 PCR Plates

- Ultra-pure polypropylene for low contamination.
- Available in 96-well or 384-well formats.
- Alphanumeric labeling for identification.
- Fits most thermal cyclers.
- Sterile, RNase- and DNase-free options.
- Skirted, semi-skirted, and non-skirted versions available.



Item	Volume/Specification	Pack Size	Catalog No.
96-Well PCR PlateNo Skirt	0.1ml, Transparent or Clear	10	CS-96-010-PCR-NS-C
96-Well PCR PlateNo Skirt	0.1ml, White	10	CS-96-010-PCR-NS-W
96-Well PCR PlateNo Skirt	0.2ml, Transparent or Clear	10	CS-96-010-PCR-NS-C
96-Well PCR PlateNo Skirt	0.2ml, White	10	CS-96-010-PCR-NS-W
96-Well PCR PlateSemi Skirt	0.2ml, Transparent or Clear	10	CS-96-020-PCR-SS-C
96-Well PCR PlateSemi Skirt	0.2ml, Black	10	CS-96-020-PCR-SS-B
96-Well PCR PlateSemi Skirt	0.2ml, White	10	CS-96-020-PCR-SS-W
96-Well PCR PlateFull Skirt	0.1ml, Transparent or Clear	10	CS-96-020-PCR-FS-C

# 8 Optical Sealing Film

- High optical clarity for real-time monitoring.
- Secure seal prevents evaporation and contamination.
- Temperature-resistant to maintain seal integrity.
- Pierceable for automated liquid handling.
- DNase/RNase-free to protect nucleic acids.




Item	Pack Size	Catalog No.
Optical Sealing Film	100	CS-PCR-OSF-100

# 9

## Deep well Plates

- Made of polypropylene
- Dnase/Rnase free
- Compatible with Automated system
- Sealing film
- Autoclavable
- Suitable for multi channel pipettes




Item	Volume	Pack Size	Catalog No.
96-Deep-Well Plate Square Well Round Conical Bottom (Elution Plate)	200ul	50	CS-96-0.5-RC
96-Deep-Well Plate Square Well Round Conical Bottom (Elution Plate)	1.6ml	50	CS-96-1.6-RC
96-Deep-Well Plate Square Well Round Conical Bottom (Elution Plate)	2.2ml	50	CS-96-2.2-SC
96-Deep-Well Plate Square Well Round U- Bottom (Elution Plate)	2.2ml	50	CS-96-2.2-SU

# 10

## Microtiter Plates

- Made of polystyrene and available in sterile or non – sterile options
- Available in various well counts including 6, 24, 96 or 384 wells.
- Different color options – clear, white or black, depending on their intended use.
- Designed for easy stacking and storage.
- Available with lids to prevent contamination
- Different well surface treatments are available for cell growth or non-stick surfaces.
- Barcode option also available for easy identification and tracking.
- Compatible with automatic pipetting systems, plate readers, and other laboratory instruments.



Item	Pack Size	Catalog No.
Microtitration Plates, Flat Bottom	100	CS-MTP-96F-S
Microtitration Plates, Round Bottom	100	CS-MTP-96U-S
Microtitration Plates with Lids, Flat Bottom	100	CS-MTP-96F-L-S
Microtitration Plates with Lids, Round Bottom	100	CS-MTP-96U-L-S



# 11

## PCR Tube Racks

- Made from high-quality materials
- Durable and built to withstand the rigors of daily laboratory use
- Available in various sizes and configurations to accommodate different types and numbers of tubes.
- Provides convenience and flexibility for laboratory professionals




Item	No. of Holes	Catalog No.
PCR Tube Rack	96	CS-PCR-R96

# 12

## Microcentrifuge Tubes Racks

- Made of autoclavable polypropylene, available in rainbow packs (B, O, R, G, Y) or single color.
- Holds multiple microcentrifuge tubes securely.
- Prevents spills and contamination.
- Compact design saves lab bench space.
- Easy-to-read labels for sample organization.
- Washable and autoclavable.
- Suitable for various lab tests and research.



Item	No. of Holes	Catalog No.
Microcentrifuge Tube Rack	96	CS-MCTR-96

# 13

## CRYO BOX

- Made from virgin polypropylene with 99.99% purity, suitable for use in IP grade applications.
- With a lid featuring alphanumeric labeling, it can be used from -90°C to 121°C.
- Autoclavable
- Available in various sizes and configurations to fit different numbers of tubes (81, 96 or 100)

Item	No. of Holes	Catalog No.
Cryobox - 81	81	CS-CB-81
Cryobox - 100	100	CS-CB-100

## Contact Us

For inquiries or to place an order, reach out to us today!

Email: [cytoscientific@gmail.com](mailto:cytoscientific@gmail.com)

Phone: [+91-8130681839](tel:+91-8130681839), [+91-7696220884](tel:+91-7696220884)

Website: [www.cytoscientific.com](http://www.cytoscientific.com)